




SPARK WAVE® THERAPY

EFFECTIVELY TREATING A WIDE RANGE
OF VETERINARY INDICATIONS

ORTHOPAEDICS / TRAUMATOLOGY / WOUND HEALING

vetgold100®



vetgold100[®] – MTS SPARK WAVE[®] TECHNOLOGY

vetgold100, MTS's premium and versatile therapy device for the treatment of a broad spectrum of veterinary indications. Intuitive operating technology, smart design and highly advanced MTS Spark Wave Technology guarantee shock wave therapy at its best.

Advanced Spark Wave Technology guarantees the most effective shock waves, consistently providing superior treatment quality and optimal treatment results.

An ergonomic touch screen and touch wheel control panel, together with a selection of lightweight handheld applicators, ensure an easy operating experience for the user.

A compact design enables **vetgold100** for portable use when mobility is required. Specially fabricated hardcases for safe and secure transfer of the device are available.



SPARK WAVE® THERAPY FOR VETERINARY INDICATIONS

Tissue regeneration with Spark Wave Therapy

Spark Waves promote sustainable regeneration of blood vessels and skin as well as of musculoskeletal and nervous tissues, which has been widely proven by comprehensive research and numerous clinical trials.

Spark Wave Therapy has a highly beneficial effect on most every kind of tissue. Cells that are treated with Spark Waves upregulate the expression of VEGF (vascular endothelial growth factor) and other growth factors, which activate neovascularisation, tissue remodelling and regeneration.

Additionally, the anti-inflammatory as well as the antibacterial effect of Spark Waves, in combination with various other biological regeneration processes, support the natural self-healing of the diseased tissue.

The efficacy of Spark Wave Therapy has been thoroughly proven in clinical studies and state-of-the-art research has demonstrated its curative effects on veterinary indications.

Treatment with **vetgold100** stimulates healing processes at the cellular level:

- Activating endogenous growth factors (e. g. BDNF, BMP, FGF-2, IGF-1, NO, PCNA, PLGF, TGF-β, VEGFs)
- Inducing stem cell proliferation and migration
- Stimulating neovascularisation and angiogenesis
- Enhancing blood circulation and cell metabolism
- Reducing cell apoptosis and tissue necrosis
- Acting anti-inflammatory and antibacterial
- Stimulating the catabolism and the natural anabolic growth function of tissues (blood and lymphatic vessels, skin, connective tissue, bones, cartilage, tendons, joints, smooth and striated muscles, nerves ...)
- Promoting tissue remodelling and regeneration

Practical selection of applicators for optimal treatment results

vetgold100 utilizes three specially designed applicators, each to treat a specific set of veterinary indications. All MTS applicators feature MTS's highly sophisticated SmarTrode® technology, which provides consistent Spark Wave performance quality.

Applicators are easily interchangeable with the device and offer the user an economical and practice-oriented treatment.

The applicator OE050 and the OE035 are perfectly suited for musculoskeletal applications. The OP155 applicator is especially effective for dermatological indications.

Unfocused applicator OP155 Soft-focused applicator OE035 Focused applicator OE050

Perfectly suited for the treatment of domestic animal and equine indications

- > Soft tissue indications
- > Tendon and ligament injuries
- > Tendinopathies, acute and chronic, e. g. tendinosis calcarea
- > Degenerative, chronic joint diseases / arthroses, e. g. hip or elbow dysplasia
- > Sesamoiditis
- > Spondylosis
- > Sacroiliac and back problems (e. g. kissing spines syndrome)
- > Problems of neck ligaments (ligamentum nuchae)
- > Pseudarthroses and delayed unions
- > Stress and avulsion fractures
- > Periostitis
- > Bone cysts (MC3)
- > Navicular syndrome
- > Interosseus ligament tears
- > Lick granuloma
- > Non-healing wounds



Benefits of vetgold100® at a glance

- Safe and effective therapy system for regenerative medicine
- Effective treatment of acute or chronic tendinopathies and musculoskeletal disorders
- Increases muscle resilience and performance capacity
- Promotes wound healing
- Non-invasive, tissue-conserving, painless and highly efficient method to enhance musculoskeletal tissue remodelling and regeneration
- Activates the metabolism and proliferation of progenitor cells
- Improves blood circulation and tissue perfusion
- Regenerative and anabolic action
- Anti-inflammatory and antibacterial effect
- Analgesic power
- User-friendly and patient-convenient
- Premium MTS Spark Wave® Technology



vetgold100®

Applicator OP155 unfocused

Max. Penetration Depth (- 6 dB):	61 mm
Max. Penetration Depth (5 MPa):	119 mm
Energy Flux Density (EFD):	0.005 – 0.176 mJ / mm ²
Acoustic Pulse Energy (E _{total} / - 6 dB):	0.40 – 3.88 mJ
Acoustic Pulse Energy (E _{total} / 5 MPa):	0.40 – 17 mJ

Applicator OE035 soft-focused

Max. Penetration Depth (- 6 dB):	60 mm
Max. Penetration Depth (5 MPa):	86 mm
Energy Flux Density (EFD):	0.005 – 0.135 mJ / mm ²
Acoustic Pulse Energy (E _{total} / - 6 dB):	0.34 – 4.22 mJ
Acoustic Pulse Energy (E _{total} / 5 MPa):	0.35 – 12.10 mJ

Applicator OE050 focused

Max. Penetration Depth (- 6 dB):	49 mm
Max. Penetration Depth (5 MPa):	60 mm
Energy Flux Density (EFD):	0.03 – 0.27 mJ / mm ²
Acoustic Pulse Energy (E _{total} / - 6 dB):	0.30 – 6.15 mJ
Acoustic Pulse Energy (E _{total} / 5 MPa):	0.61 – 19.50 mJ

TECHNICAL DATA

Dimensions:	21.8 x 40 x 45.9 cm
Weight:	18 kg
Power Supply:	100 - 240 VAC, 50 - 60 Hz
Power Input:	200 VAC
Frequency:	0.5 - 8 Hz

MTS Medical specialises in the research, development and manufacturing of Spark Wave® Therapy systems. Our therapy systems are successfully used worldwide in the healthcare fields of orthopaedics, physiotherapy, traumatology, urology, wound care and aesthetics. MTS therapy systems utilize MTS's own patented Spark Wave® Technology, which guarantees the most effective shock waves for pain-free, high quality treatment as well as quick and excellent treatment results.

Let yourself be convinced by the Spark of Evidence.



ESWT



MADE IN



GERMANY

MTS Medical UG
Robert-Bosch-Str. 18
D-78467 Konstanz
+49 (0) 7531 36185-0
sales@mts-medical.com
www.mts-medical.com
www.mts-science.com



MTS

SPARK OF EVIDENCE